


ACE-801
Front closure assembly:
chafe, bar, two screws



ACE-802
Hinged pubic assembly:
bar, plate, pad, padding,
screws




ACE-803
Hinged sternal
assembly: bar, plate,
pad, padding, screws




ACE-804
PV assembly: bar,
plate, Velcro, pair of
PV pads (ACE-806),
screws



ACE-805
PS assembly: bar, V
slotted bar, pair of PS
pads (ACE-807),
screws (screws not
shown)



ACE-806
PV pad assembly: pad,
padding, Neoprene disc,
Velcro coin (pair of 2)



ACE-807
PS pad assembly: pad,
padding, silicone hinge,
retainer plates, screws
(pair of 2) (retainer
plates and screws not
shown)



ACE-808
Side pad assembly: bar,
plate, padding, screws
(screws not shown)



ACE-809
Back strap
assembly: strap,
Y-tab, back pad,
padding



ACE-810
Back strap
assembly, large:
strap, Y-tab, large
back pad, padding



ACE-811
Back Pad, regular:
polyethylene plate,
padding (8.25" x 4.5")




ACE-812
Back Pad, large:
polyethylene plate,
padding (10" x 5.5")




ACE-813
Back strap: long
strap, Y-tab
(universal size)



ACE-948
Y-tab (part of
ACE-813)




ACE-814
Foam padding,
complete set: for
ACE-ST



ACE-814
Foam padding,
complete set: for ACE-
PV and ACE-PS



ACE-910
Vertical bar,
regular



ACE-912
Vertical bar, long (adds 3.25"
to height of ACE Brace)



ACE-911
Horizontal bar,
regular



ACE-913
Horizontal bar,
long (adds 1.5"
to width of ACE
Brace)