

NINJA™ Ordering Information

NINJA
NINJA Belt

Belt Only



NJLSM
NINJA Mid LSO

Adjustable
Lordosis



NJLSO
NINJA LSO

Adjustable
Lordosis &
Kyphosis






NJPRO
NINJA PRO


Molded
Lordosis
05°, 15°, 25°



NINJA Heights



NINJA Size Chart





NINJA Front Heights

XX - 5X (8")

SM - XL (7 1/2")

XX - 5X (6 1/2")

SM - XL (6")



NINJA Ordering Grid				
Type	Belt	Size	Front	Back
NINJA	LP	SM	MPE	05°
NJLSM	SP	MD	KYD	15°
NJLSO		LG	PKYD	25°
NJPRO		XL		
		XX		OSR
		3X		
		4X		
		5X		

N J P R O - S P - L G - M P E - 1 5

Example: NINJA PRO, Standard Belt, Large, Polyethylene Insert , with 15° Lordosis.

Open Spinal Relief



Open Spinal Relief provides an oval spinal cutout over the spine, allowing visual wound inspection.

Accessories

NINJA Anterior Inserts



MPE -
Polyethylene Insert
(Standard)
Not Contoured
No Extender Holes



KYD -
Molded Kydex
Extender Ready
PKYD -
Pendulous KYD

The Extender™



EXT STP
White Straps
Soft Cover



EXT STP
No Straps
No Soft Cover

Mickey



Two posterior thoracic options provide over-the-shoulder strapping alternatives. Both attach to the back of a spinal orthosis by Velcro.

The **Mickey™** contains its own straps that start posteriorly, pass over the shoulders then back under the axilla, and reattach to itself.

Minnie



The **Minnie™** has no straps, but receives the Extender's straps, which pass either over the shoulders or under the axillae and attach to the Minnie by Velcro.

The Extender™

The Extender may be attached to a NINJA with the optional Extender ready anterior plate. There are many Extender variations:

1. Six styles of sternal/pectoral pressure plates
2. Fixed or Adjustable Extension
 - Both with sliding height adjustments
3. Black, white or no shoulder straps
4. With or without trilaminar soft pad covers

For more Extender options, please see our Extender Guide