

# OrthoLux™ **Ordering Information**

OrthoLux Ordering Grid				
OrthoLux style	sex	size	front	back
		SM		
		MD	NT	
		LG	PA	00
		XL	MP	07
LSM		XX	PM	15
LSO		3X	LP	25
TLS	U	4X	PL	35

**ORT - LSO - U - LG - NT - 15**

Example: LSO, unisex, large, neutral front, 15° lordosis back

OrthoLux Unisex Size Chart			
SM	Hip 29" - 35"	Waist 28" - 34"	Rib 28" - 34"
MD	Hip 32" - 39"	Waist 30" - 37"	Rib 30" - 36"
LG	Hip 36" - 42"	Waist 34" - 42"	Rib 34" - 40"
XL	Hip 39" - 47"	Waist 38" - 45"	Rib 38" - 45"
XX	Hip 45" - 52"	Waist 43" - 51"	Rib 43" - 51"
3X	Hip 51" - 56"	Waist 47" - 55"	Rib 47" - 55"
4X	Hip 57" - 63"	Waist 53" - 61"	Rib 53" - 61"

OrthoLux Front Heights			
	Low-Profile (LP/PL)	Mid-Profile (MP/PM)	Neutral (NT/PA)
SM	LP/PL* 6 1/2"	MP/PM* 7 1/2"	NT/PA* 9 3/8"
MD	LP/PL* 7"	MP/PM* 8 1/2"	NT/PA* 10 1/4"
LG	LP/PL* 7"	MP/PM* 8 1/2"	NT/PA* 10 3/8"
XL	LP/PL* 7 3/8"	MP/PM* 9"	NT/PA* 11 1/4"
XX	LP/PL* 7 3/8"	MP/PM* 9 3/8"	NT/PA* 11 1/4"
3X	LP/PL* 7 3/4"	MP/PM* 9 3/8"	NT/PA* 11 1/4"
4X	LP/PL* 7 3/4"	MP/PM* 9 3/8"	NT/PA* 11 1/4"

*(These are fabric heights; insert heights are about 2" less on average)*

\*PL = Pendulous Low Profile; PM = Pendulous Mid Profile; PA = Pendulous Xiphoid Height

OrthoLux Back Heights			
	LSO Mid-Profile (LSM)	LSO (LSO)	TLSO (TLS)
SM	LSM 11 1/2"	LSO 15 1/4"	TLSO 18 1/4"
MD	LSM 11 3/4"	LSO 15 3/4"	TLSO 18 3/4"
LG	LSM 12 1/4"	LSO 16"	TLSO 19"
XL	LSM 12 1/4"	LSO 16"	TLSO 19"
XX	LSM 13 1/8"	LSO 18"	TLSO 21"
3X	LSM 13 1/8"	LSO 18"	TLSO 21"
4X	LSM 13 1/8"	LSO 18"	TLSO 21"

*(These are fabric heights; insert heights are about 2" less on average)*

The OrthoLux Orthosis is available in three front heights (each with 2 contours - neutral or pendulous), three back heights (each with five lordoses - 00°, 07°, 15°, 25° and 35°), and seven sizes. The OrthoLux is a completely modular system, so you can mix different sizes or styles of back and front for maximum fitting precision. (See charts for exact sizes and heights.)

## ACCESSORIES

The "Mickey" is a TLSO add-on that applies a moderate extensive force to the thoracic spine. The "Mickey" can be Velcroed to the upper back panel of the OrthoLux. The straps pass over the shoulder, then extend downward anterior to the pectoral muscle and then pass back through the axillae to attach posteriorly.



The lateral extension panels can be used to either fill the lateral gap between the front and back panels, or to add up to 8" to the circumference of the OrthoLux.

## EXTENDERS



Extender EXT-STP-ADJ-WHT-WSC

Pictured above is our most popular Extender™, the EXT-STP (shown with adjustable extension, over-the-shoulder straps, and soft pad cover). There are six variations of Extender tops, and each is available with three options: Adjustable Extension, White or Black Straps, and Soft Pad Cover(s). (For detailed Extender ordering information, please see the Extender ordering Guide.)

Copyright 2007  
HOPE Orthopedic  
ph: 800-613-8852  
fx: 888-440-1217

