

OrthoLux Lite™ **Ordering Information**

OrthoLux Lite Ordering Grid

OrthoLux style	sex	size	front	back
		SM	NT	
		MD	PA	00
		LG	MP	07
LSM		XL	PM	15
LSO		XX	LP	25
TLS	U	3X	PL	35

OLL - LSO - U - LG - NT - 15

Example: LSO, unisex, large, neutral front, 15° lordosis back

OrthoLux Lite Unisex Size Chart

SM	Waist Circumference	26" - 32"
MD	Waist Circumference	30" - 37"
LG	Waist Circumference	34" - 41"
XL	Waist Circumference	38" - 47"
XX	Waist Circumference	44" - 54"
3X	Waist Circumference	50" - 62"

Although the size ranges to the left specify "Waist Circumference," since this is an all-elastic orthosis, it is best to measure the *largest* circumference, whether that be at the level of the waist or of the hip.

OrthoLux Lite Front Heights

	Low-Profile (LP/PL)	Mid-Profile (MP/PM)	Neutral (NT/PA)
SM	LP/PL* 6"	MP/PM* 8"	NT/PA* 8 3/4"
MD	LP/PL* 6"	MP/PM* 8"	NT/PA* 8 3/4"
LG	LP/PL* 6"	MP/PM* 8"	NT/PA* 8 3/4"
XL	LP/PL* 6 1/2"	MP/PM* 8 1/2"	NT/PA* 9 1/4"
XX	LP/PL* 6 1/2"	MP/PM* 8 1/2"	NT/PA* 9 1/4"
3X	LP/PL* 6 1/2"	MP/PM* 8 1/2"	NT/PA* 9 1/4"

(These are fabric heights; insert heights are about 2" less on average)

*PL = Pendulous Low Profile; PM = Pendulous Mid Profile; PA = Pendulous Xiphoid Height

OrthoLux Lite Back Heights


	LSO Mid-Profile (LSM)	LSO (LSO)	TLSO (TLS)
SM	LSM 13 1/4"	LSO 15 3/8"	TLSO 18 3/8"
MD	LSM 13 1/4"	LSO 15 3/8"	TLSO 18 3/8"
LG	LSM 13 1/4"	LSO 15 3/8"	TLSO 18 3/8"
XL	LSM 13 3/4"	LSO 15 3/4"	TLSO 18 3/4"
XX	LSM 13 3/4"	LSO 15 3/4"	TLSO 18 3/4"
3X	LSM 13 3/4"	LSO 15 3/4"	TLSO 18 3/4"

(These are fabric heights; insert heights are about 2" less on average)

The OrthoLux Lite Orthosis is available in three front heights (each with 2 contours - neutral or pendulous), three back heights (each with four lordoses - 00°, 07°, 15°, and 25°), and six sizes. The OrthoLux is a modular system, so you can mix different sizes or styles of back and front for maximum fitting precision. (See charts for exact sizes and heights.)

ACCESSORIES

The "Mickey" is a TLSO add-on that applies a moderate extensive force to the thoracic spine. The "Mickey" can be Velcroed to the upper back panel of the OrthoLux. The straps pass over the shoulder, then extend downward anterior to the pectoral muscle and then pass back through the axillae to attach posteriorly.



EXTENDERS



Pictured above is our most popular Extender™, the EXT-STP (shown with adjustable extension, over-the-shoulder straps, and soft pad cover). There are six variations of Extender tops, and each is available with three options: Adjustable Extension, White or Black Straps, and Soft Pad Cover(s). (For detailed Extender ordering information, please see the Extender ordering Guide.)

Copyright 2007
HOPE Orthopedic
ph: 800-613-8852
fx: 888-440-1217

HOPE
ORTHOPEDIC
www.hopeortho.com