

TorsoLux™ *Ordering Information*

TorsoLux Ordering Grid				
TorsoLux	style	size	front	back
		S M		
		M D		
		LG		
		X L	SP	
		X X	P A	0 5
		3 X	M P	1 5
	LSO	4 X	P M	2 5
T S L	- L S O	- L G	- S P	- 1 5

Example: LSO, large, standard profile front, 15° lordosis back

TorsoLux Unisex Size Chart			
SM	Hip 29" - 35"	Waist 28" - 34"	Rib 28" - 34"
MD	Hip 32" - 39"	Waist 30" - 37"	Rib 30" - 36"
LG	Hip 36" - 42"	Waist 34" - 42"	Rib 34" - 40"
XL	Hip 39" - 47"	Waist 38" - 45"	Rib 38" - 45"
XX	Hip 45" - 52"	Waist 43" - 51"	Rib 43" - 51"
3X	Hip 51" - 56"	Waist 47" - 55"	Rib 47" - 55"
4X	Hip 57" - 63"	Waist 53" - 61"	Rib 53" - 61"

TorsoLux Front Heights			
	Mid-Profile (MP/PM)		Neutral (SP/PA)
SM	MP/PM* 7"		SP/PA* 9"
MD	MP/PM* 7"		SP/PA* 9"
LG	MP/PM* 7"		SP/PA* 9"
XL	MP/PM* 7"		SP/PA* 9"
XX	MP/PM* 7"		SP/PA* 9"
3X	MP/PM* 7"		SP/PA* 9"
4X	MP/PM* 7"		SP/PA* 9"

(These are fabric heights; insert heights are about 1.5" less on average)

*PM = Pendulous Mid Profile; PA = Pendulous Standard Profile (Xiphoid Height)

TorsoLux Back Heights		
	LSO Plastic Frame	LSO Foam Padding
SM	LSO 14"	LSO 15 1/4"
MD	LSO 14"	LSO 15 1/4"
LG	LSO 14"	LSO 15 1/4"
XL	LSO 14"	LSO 15 1/4"
XX	LSO 14"	LSO 15 1/4"
3X	LSO 14"	LSO 15 1/4"
4X	LSO 14"	LSO 15 1/4"

The TorsoLux Orthosis is available in two front heights (each with 2 contours - neutral or pendulous), one back height* (LSO) with three lordoses (05°, 15°, and 25°), and seven sizes. The TorsoLux is a modular system, so you can mix different sizes of back and with different styles and sizes of front for maximum fitting precision. (See charts for exact sizes and heights.) *TLSO and/or LSM posterior heights are under consideration. Please let us know if they would be helpful to you.

ACCESSORIES

Because some orthotists prefer to stock modular orthoses by components and then "build" them to order, the TorsoLux Orthosis is available in components. You may order the front fabric and/or panel separately, but keep in mind that the back frame assembly will always come with the straps, as they are riveted to the back polymer frame. However, the posterior tri-laminate padding is always available separately.

The anterior polymer panels (standard profile & mid profile) are pre-drilled to accept any model of our Pectoral or Sternal Extenders. Thus the orthotist can easily make a last-minute decision to add an Extender while at the hospital, and can attach, apply and adjust the Extender with the tools included.

EXTENDERS



Extender EXT-STP-ADJ-WHT-WSC

Pictured above is our most popular Extender™, the EXT-STP (shown with adjustable extension, over-the-shoulder straps, and soft pad cover). There are six variations of Extender tops, and each is available with three options: Adjustable Extension, White or Black Straps, and Soft Pad Cover(s). (For detailed Extender ordering information, please see the Extender ordering Guide.)

Copyright 2007
HOPE Orthopedic
ph: 800-613-8852
fx: 888-440-1217

